

---

# Bölüm 33

## Özefagiyal Manometre ve Mekanik Ventilasyon Sırasında Yatakbaşı Görüntüleme

---

Çeviri: Dr. Yağmur Kara, Dr.Hülya Sungurtekin

- **Giriş**
- **Özefagus Basıncı**
  - Basıncın Plevral Boşluğa İletimi
  - PEEP Titrasyonu
  - Hasta vs Ventilatör Solunum İşİ
  - Spontan Solunum ile Oto-PEEP
- **İntra-Abdominal Basınç**
  - Gastrik Basınç
  - Mesane Basıncı
- **Akciğer Volümü**
- **Yatakbaşı Görüntüleme**
  - Ultrasonografi
  - Elektriksel İmpedans Tomografi
- **Hatırlanacak Noktalar**
- **Ek Okumalar**

**Hatırlanacak Noktalar**

- Özefagus basıncı, plevral basınç yerine kullanılabilir.
- Özefagus basıncı, PEEP ayarını, hastanın solunum işini, spontan solunum sırasındaki oto-PEEP'i ve göğüs duvarı kompliyansını değerlendirmek için kullanılabilir.
- İntra-abdominal basıncın ölçülmesi, transdiyafragmatik basıncı ve intra-abdominal basıncın göğüs duvarı kompliyansına etkisini değerlendirmede faydalı olabilir.
- EELV, azot washout tekniği kullanılarak mekanik ventilasyon sırasında ölçülebilir.
- Ultrasonografi akciğer parankimini ve diyafragma fonksiyonunu değerlendirmek için kullanılabilir.
- EIT, PEEP titrasyonu sırasında rekrutment ve overdistansiyonu değerlendirmek için kullanılabilir.

**Ek Okumalar**

- Akoumianaki E, Maggiore SM, Valenza F, et al.** The application of esophageal pressure measurement in patients with respiratory failure. *Am J Respir Crit Care Med.* 2014;189(5):520-531. **Bellani G, Mauri T, Pesenti A.** Imaging in acute lung injury and acute respiratory distress syndrome. *Curr Opin Crit Care.* 2012;18(1):29-34. bu ayrı bir referans oldu için diğer referanslar gibi satırbaşı yapılarak yazılacak
- Bellani G, Rouby JJ, Constantin JM, Pesenti A.** Looking closer at acute respiratory distress syndrome: the role of advanced imaging techniques. *Curr Opin Crit Care.* 2017;23(1):30-37.
- Bikker IG, Preis C, Egal M, Bakker J, Gommers D.** Electrical impedance tomography measured at two thoracic levels can visualize the ventilation distribution changes at the bedside during a decremental positive end-expiratory lung pressure trial. *Crit Care.* 2011;15(4):R193.
- Blankman P, Hasan D, Erik G, Gommers D.** Detection of 'best' positive end-expiratory pressure derived from electrical impedance tomography parameters during a decremental positive end-expiratory pressure trial. *Crit Care.* 2014;18(3):R95.
- Bouhemad B, Brisson H, Le-Guen M, Arbelot C, Lu Q, Rouby JJ.** Bedside ultrasound assessment of positive end-expiratory pressure-induced lung recruitment. *Am J Respir Crit Care Med.* 2011;183(3):341-347.
- Bouhemad B, Zhang M, Lu Q, Rouby JJ.** Clinical review: bedside lung ultrasound in critical care practice. *Crit Care.* 2007;11(1):205.
- Branson RD, Johannigman JA.** Innovations in mechanical ventilation. *Respir Care.* 2009;54(7):933-947.
- Dellamonica J, Lerolle N, Sargentini C, et al.** PEEP-induced changes in lung volume in acute respiratory distress syndrome. Two methods to estimate alveolar recruitment. *Intensive Care Med.* 2011;37(10):1595-1604.
- Doerschug KC, Schmidt GA.** Intensive care ultrasound: III. Lung and pleural ultrasound for the intensivist. *Ann Am Thorac Soc.* 2013;10(6):708-712.
- Fumagalli J, Berra L, Zhang C, et al.** Transpulmonary pressure describes lung morphology during decremental positive end-expiratory pressure trials in obesity. *Crit Care Med.* 2017;45(8):1374-1381.

- Gardelli G, Feletti F, Nanni A, Mughetti M, Piraccini A, Zompatori M.** Chest ultrasonography in the ICU. *Respir Care.* 2012;57(5):773-781.
- Gattinoni L, Chiumello D, Carlesso E, Valenza F.** Bench-to-bedside review: chest wall elastance in acute lung injury/acute respiratory distress syndrome patients. *Crit Care.* 2004;8(5):350-355.
- Grasso S, Terragni P, Birocco A, et al.** ECMO criteria for influenza A (H1N1)-associated ARDS: role of transpulmonary pressure. *Intensive Care Med.* 2012;38(3):395-403.
- Gulati G, Novero A, Loring SH, Talmor D.** Pleural pressure and optimal positive end-expiratory pressure based on esophageal pressure versus chest wall elastance: incompatible results\*. *Crit Care Med.* 2013;41(8):1951-1957.
- Hess DR.** Approaches to conventional mechanical ventilation of the patient with acute respiratory distress syndrome. *Respir Care.* 2011;56(10):1555-1572.
- Hess DR.** Respiratory mechanics in mechanically ventilated patients. *Respir Care.* 2014;59(11):1773-1794.
- Hess DR, Bigatello LM.** The chest wall in acute lung injury/acute respiratory distress syndrome. *Curr Opin Crit Care.* 2008;14(1):94-102.
- Kobylianskii J, Murray A, Brace D, Goligher E, Fan E.** Electrical impedance tomography in adult patients undergoing mechanical ventilation: a systematic review. *J Crit Care.* 2016;35:33-50.
- Lecamwasam HS, Hess D, Brown R, Kwolek CJ, Bigatello LM.** Diaphragmatic paralysis after endovascular stent grafting of a thoracoabdominal aortic aneurysm. *Anesthesiology.* 2005;102(3):690-692.
- Loring SH, O'Donnell CR, Behazin N, et al.** Esophageal pressures in acute lung injury: do they represent artifact or useful information about transpulmonary pressure, chest wall mechanics, and lung stress? *J Appl Physiol* (1985). 2010;108(3):515-522.
- Matamis D, Soilemezi E, Tsagourias M, et al.** Sonographic evaluation of the diaphragm in critically ill patients. Technique and clinical applications. *Intensive Care Med.* 2013;39(5):801-810.
- Mauri T, Yoshida T, Bellani G, et al.** Esophageal and transpulmonary pressure in the clinical setting: meaning, usefulness and perspectives. *Intensive Care Med.* 2016;42(9):1360-1373.
- Pesenti A, Musch G, Lichtenstein D, et al.** Imaging in acute respiratory distress syndrome. *Intensive Care Med.* 2016;42(5):686-698.
- Pirrone M, Fisher D, Chipman D, et al.** Recruitment maneuvers and positive end-expiratory pressure titration in morbidly obese ICU patients. *Crit Care Med.* 2016;44(2):300-307.
- Rouby JJ, Arbelot C, Brisson H, Lu Q, Bouhemad B.** Measurement of alveolar recruitment at the bedside: the beginning of a new era in respiratory monitoring? *Respir Care.* 2013;58(3):539-542.
- Talmor D, Sarge T, Malhotra A, et al.** Mechanical ventilation guided by esophageal pressure in acute lung injury. *N Engl J Med.* 2008;359(20):2095-2104.
- Walsh BK, Smallwood CD.** Electrical impedance tomography during mechanical ventilation. *Respir Care.* 2016;61(10):1417-1424